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@neubitam

BİTAM

A photograph of a modern, multi-story building with a glass facade, situated on a hillside. The building has a prominent sign on its upper section that reads "BİTAM" in blue capital letters. The building's facade is composed of a grid of glass panels, reflecting the surrounding environment. The building is set against a backdrop of a hillside covered in sparse vegetation and trees. The sky is filled with soft, colorful clouds, suggesting a sunset or sunrise. The overall scene is captured in a wide-angle shot, emphasizing the building's scale and its integration with the natural landscape.

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SCIENCE TECHNOLOGY
RESEARCH AND APPLICATION CENTER

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Dear Researchers,

Science, the common heritage of humanity, has been shaped by the ideas produced by human beings with the effort of recognizing oneself and the universe, and the hard work done under the influence of the conditions developed by these ideas. This process will build up the knowledge step by step and transferred to today's world, it will undoubtedly continue to be the light of humanity in the future too. It is of great importance for the future of all countries in the world, especially in our country, to support the purposeful, conscious, and responsible works carried out with the aim of being beneficial to humanity and to provide the necessary conditions for their fulfillment.

Science and Technology Research and Application Center (BITAM) is dedicated to our universities, public institutions, and industry, and I believe that it will open new horizons by increasing the quality of the accumulation of knowledge to be passed on to future generations.

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BİTAM

Dear Researchers,

Necmettin Erbakan University Science and Technology Research and Application Center (BITAM), which aims to be a training center for your knowledge that you will transfer to future generations within the framework of universal values, is equipped with up-to-date and high-tech devices. Numerous analysis required for the emergence of technological, unique and national products with national resources can be performed in our unit. In addition, our laboratory, which includes many thematic area laboratories determined in line with the needs of our country in recent years, serves in a closed area of 3000 m².

We believe that our center will contribute to our valuable researchers and producers who have taken it upon themselves to build a better future for our nation. We believe that we will succeed with you in our unit where researchers from different disciplines are together.

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
SURFACE AND DIMENSIONAL CHARACTERIZATION INSTRUMENTS

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84-85

3D TASARIM VE PROTOTİPLEME CİHAZLARI

3D DESIGN AND PROTOTYPE INSTRUMENTS

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Device Model : Thermo Scientific / Nicolet IS20

With advanced LightDrive optical engine technology, the device identifies pollutants in the environment or characterizes small amounts of errors with its high spectral resolution ($0,25 \text{ cm}^{-1}$) and single noise ratio (50,000: 1). It has a better spectral resolution than $0,25 \text{ cm}^{-1}$.

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Infrared (IR) Spectroscopy is a spectroscopy method based on the absorption of infrared light by matter and obtaining the emission infrared spectrum. It provides the recognition of the functional groups in the structure by measuring the vibration frequencies of the bonds in the molecules.

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In the structural analysis of carbohydrates, proteins and amino acids,

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Structure determination and characterization of compounds such as plastics, fillers, paints, rubbers, coatings, resins and adhesives.

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In the determination of molecular interactions.

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Device Model : Agilent 7900 ICP-MS

Thanks to its ability to analyze many elements, at the same time, it is widely used in qualitative analysis and determination of isotope ratios, as well as quantitative and semi - qualitative determinations of the majority of the elements in the periodic table, especially metallic elements. Trace element analysis can be done at the ppb level.

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Inductively coupled plasma-mass spectrometry is an analysis method in which the molecular bonds are broken and the atoms are ionized by sending the sample to a high temperature plasma. It usually happens by ionizing the sample in solution with inductively coupled plasma (Ar plasma) and then determining the mass / charge ratios using a mass spectrometer to separate and measure these ions.

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One of the most common uses is environmental analysis (Soil samples, natural spring water, drinking water, wastewater, etc.). It is also frequently used for metal and trace element analysis in the fields of mining, geology, petrochemicals, food, drug and toxicology.

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The determination limit of many elements for ICP-MS is less than one billionth (ppb). Therefore, it is the most suitable method for trace element determination and quantification.

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Since the solution ICP-MS analysis technique is applied in the analyzes, solid samples are converted into solutions before analysis. Different sample preparation methods are applied according to the type of the sample (geological, environmental, biological etc.) and the elements to be analyzed.

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Device Model : GNR S7 METAL LAB Plus

Optical emission spectrometry (OES) is an analysis technique used especially in qualitative and quantitative characterization of metals, alloys. OES can analyze a wide range of elements from Lithium to Uranium in solid metal examples covering a wide concentration range, giving very high accuracy, high precision and low detection limits.

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This method is based on the principle that each element in the molecule emits characteristic spectral lines at certain frequencies as a result heating the sample with flame. In order to distinguish the lines from each other, a grid reflector is used, and their intensity is measured by a sensor. Thus, the elements in the sample and the percentage of amounts of these elements are determined.

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Alloys containing Iron (Fe), Aluminum (Al), Copper (Cu), Nickel (Ni), Zinc (Zn) and Titanium (Ti).

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Computer controlled and fully automated combustion source. Computer controlled interchangeable combustion parameters for special applications.

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Due to its large optical tank with CCD detector with 500 mm focal length, it has low interference effect, high resolution, analysis stability and repeatability.

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Device Model : Setaram Labsys Evo

Thermogravimetric analysis (TGA) is an analytical technique used to determine a material's thermal stability and its fraction of volatile components by monitoring the weight change that occurs as a sample is heated at a constant rate. Energy behavior of reaction types that occur during the heating of the sample can be described by differential scanning calorimetry analysis (DSC).

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In thermogravimetric analysis, before the analysis, the deviated sample is recorded by reading the deviations that occur with heating. Differential scanning calorimetry is a thermoanalytical technique in which the difference in the amount of heat required to increase the temperature of a sample and reference, measured as a function of temperature. The sample and reference are kept at the same temperature throughout the experiment.

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Monitoring of the mass change in materials which is caused by dehydration, degradation, oxidation, and corrosion along with the thermal change.

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Differential Scanning Calorimetry

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Determination of the amount of energy change of reactions occurring with thermal change in materials.

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Thermal stability of materials, decomposition temperature and rate can be determined by heating.

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Device Model : SHIMADZU UV-3600 PLUS

It has 3 detectors consisting of PMT (photomultiplier tube) in the Ultraviolet and Visible Region, and InGaAs and cooled PbS detectors in the near infrared region. Measurements can be performed with high sensitivity over the whole measurement range due to the InGaAs detector. (0.00003 Abs noise level at 1500 nm). It can be operated over the wide range of wavelengths, i.e. 185-3300 nm.

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It is a spectroscopy method based on the principle of measuring the decrease in the intensity of a beam reflected from the sample surface or passing through the sample in the ultraviolet, visible and near infrared (UV-Vis-NIR) region. The difference in the light absorption properties of each substance is a distinctive criterion for molecules.

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It is used to identify organic molecules, ions or complexes.

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Water (drinking water, wastewater, etc.) analysis

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It is used to determine the absorption bands formed by the interaction of nanomaterials with light at different wavelengths.

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Device Model : Rigaku NEX CG

X-Ray Fluorescence Spectrometer (XRF) Sodium (Na) - Uranium (U) is used to determine the chemical composition of solid, liquid and powder samples qualitatively and quantitatively in terms of their elements and oxides. Precise analysis from ppm to percent level can be carried out with the selection of appropriate materials and programs.

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As a result of the excitation of the atoms in a material with x-rays, their electrons in close orbits reach higher energy levels. As the excited electrons fall back to their lower energy level, they give back the excess energy in the form of light. By determining these characteristic photon energies, their chemical composition is determined qualitatively and quantitatively.

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It is used in metal, concrete, polymer, plastic, mining, metallurgy, mineral, geology and food industry.

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The measurement time varies from seconds to 30 minutes, depending on the element to be measured and the number of repetitions.

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The sample is not damaged during the analysis.

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Device Model : Bruker, 400 MHz, Avance NEO (Broadband)

NMR is an extremely powerful analytical method and an ideal tool to address a large range of analytical questions.

Structure Elucidation, Verification and Molecular Interaction

- Examine functional groups and how atoms are connected within a molecule.
- Study molecules at atomic resolution

Uncover molecular interactions at the atomic level.

NMR selectively detects nuclei such as ^1H , ^{13}C , ^{19}F and ^{31}P which are present in organic molecules. This is not apparent with other analytical techniques.

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The principle behind NMR is that NMR active nuclei have spin and all nuclei are electrically charged. If an external magnetic field is applied, an energy transfer is possible between the base energy to a higher energy level (generally a single energy gap). The energy transfer takes place at a wavelength that corresponds to radio frequencies and when the spin returns to its base level. The signal that matches this transfer is measured and processed to yield an NMR spectrum. By this method, NMR measurement can be determined.

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Instrument Properties ve Application Area

1

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Analysis and Identification of Molecular Structure

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A wide range of applications including organic chemistry, inorganic chemistry, biochemistry, drug chemistry, new materials, petrochemistry, etc.

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NMR can be able to make ingredient analysis for the field of food chemistry, biochemistry, physiology and material sciences

3

Kantitatif ve Fiziksel Analiz

Quantitative and Physical Analysis

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Some physical and chemical properties of the compound can be examined in such a field of Chemistry.

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Device Model : Panalytical/ Empyrean

It is based on the principle that each crystal phase refracts X-rays in a characteristic order depending on its atomic sequences. For each crystalline phase, these diffraction profiles define that crystal like a fingerprint. The X-Ray Diffraction analysis method does not destroy the sample during analysis and enables the analysis of even a small amount (liquid, powder, crystal and thin film).

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With the applied method, the necessary peak corrections are made on the patterns taken from the high resolution Pixel based detector, using databases containing crystal structure files, and the standard deviation rates and the percent ratios of the minerals are given. The quantitative ratios of minerals are calculated using the Rietveld method.

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Determination of phases in powder and solid samples

2

İnce Film

Thin Film

İnce film faz ve Rietveld analizi

Thin film phase and Rietveld analysis

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Determination of crystallographic structure

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Device Model : Agilent Cary Eclipse

Agilent Cary Eclipse equipped with a microplate reader (96 wells) has the xenon (Xe) flash lamp enables highly sensitive, high signal-to-noise ratio and kinetic researches. The Cary Eclipse can scan at a rate of 24,000 nm/min without peak shift using Xe lamp technology. With the microplate reader, it provides wide wavelength scanning in less than 30 seconds.

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It is a spectroscopy method based on the excitation of electrons in molecules and the measurement of the characteristic wavelength of the fluorescent light emitted. It is used effectively in the identification and characterization of small amounts of substances.

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*Fluorescence, phosphorescence
and luminescence measurements.*

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*Determination of amino acids in
biological samples.*

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Surface and Size Characterization Instruments

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Device Model : Quantachrome – Quadrasorb Evo 4

By using the BET device, surface area measurements, micro, meso and macropore size and pore size distribution can be determined by physical adsorption method in solid or powder samples.

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Surface area and porosity measurement is carried out based on the nitrogen (N_2) gas adsorption technique at 77 K liquid nitrogen environment by Brunauer, Emmett and Teller (BET) method.

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It has 4 independent analysis stations used in surface area and pore size measurements.

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BET surface area and pore size distribution can be reached reliably and accurately due to the adsorption and desorption capacities.

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Device Model : Micromeritics / Autopore IV 9500

It is used to determine pore size, pore size distribution and surface area measurements and mass density in powder or bulk samples. It is used to calculate the porosity and dimensionless volume distribution of materials such as ceramics, mineral and mineral products, building materials, biological implants, membranes, pharmaceutical raw materials, metallurgical powders, abrasives and polymers.

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The measurement is based on the principle that a non-reactive and non-wetting liquid, such as mercury, cannot enter fine pores unless sufficient pressure is applied. The device determines the amount of mercury decreasing in the cell against the applied pressure from the change of capacitance between the cell ends. Pore size is calculated as a function of pressure.

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Yoğunluk

Density

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Determination of relative density

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Device Model : Micromeritics – Accupyc II 1340

It is used to determine the volume and density value of various materials (such as metal, ceramic polymer and foam) with a maximum volume of 100 cm³ using helium gas. The resolution of the sensitive digital pressure gauge is 0.0001 psi. Due to the high number of pores in which He gas will effect the material, density determination with helium pycnometer is more reliable and sensitive than the other methods.

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It aims to find volume and real density by using Archimedes' fluid overflow principle and Boyle's Law. For the accuracy of the measurements to be maximum, the overflowing fluid must be an inert gas that can enter all pores except the smallest ones. Therefore, helium gas approaching to indentations and pores with a diameter of 0.25 mm is considered appropriate due to its small atomic dimensions.

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The volume of the material is determined with the helium gas.

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The density of the material is calculated by adding the mass value.

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Device Model : Micromeritics- Nanoplus 3

Almost all the surface of the macroscopic materials has a electrical charge when they are in the form of particles or in contact with a liquid. Zeta potential is a critical and useful indicator to measure the stability of the colloidal suspensions or emulsions. The higher zeta potential makes the suspension stable, and thus, overcomes the natural clustering trend.

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Dynamic light scattering, also called photon correlation spectroscopy and electrophoretic light scattering techniques are used to determine the particle size and zeta potential. The zeta potential are determined by indicating the speed of the particles that are moved under a certain potential. The particle size is determined through measuring the intensity of the scattered light by Brownian motion.

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Polymers and proteins, pigments, wastewater treatment, ceramic processes, liquid ink and toners etc.

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To determine the size of nanoparticles in the fields of chemistry, medicine, pharmacy and materials engineering etc.

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Measures particle size of the colloidal samples in the range of 0.1 nm and 12.3 μ m

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Device Model : Cilas 1190

The Cilas particle size analyzer provides fast, accurate and easy particle distribution with both dry and wet dispersion units. It has a feature to measure a size range from nanometers to millimeters. Three lasers with extremely high accuracy and repeatability are used to measure the widest particle range. In liquid measurements, particles are dispersed by ultrasonic effect.

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Mie Theory is valid in the size analysis of small particles while the Fraunhofer theory, which does not take into account the diffusion and absorption, is valid for large particles. Mie's theory expresses the refraction and permeability of the light around the particle in the medium. The refractive indices of both the sample and the medium must known in order to use the Mie model.

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Dry or wet measurements with high sensitivity in a wide measuring range / 0.1-2500 μm in dry mode, 0.04-2500 μm in wet mode

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In the fields of ceramics, powder metallurgy, mineral processes, powdered food products, inks, etc.

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
Device Model : Biolin Scientific Attension- Theta Lite

The surface/interface tension is an optical tensiometer system that can measure the shape and volume of the droplet created with free surface energy. In optical tensiometer systems, the evaporation amount of the droplet can be minimized by using a cold monochromatic LED light source. Information about the wettability of liquids and fluid spill on the surface can be obtained by this technique.

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The contact angle device analyzes the changes in the form of the drops over time by recording the images of the drop. The drop shape is a function of the surface tension of the liquid and consequently the difference in density between the liquid and the medium in which the liquid is located. On solid surfaces, the dynamic and static contact angle varies depending on the free surface energy of the solid.

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Static contact angle, Dynamic contact angle, Free surface energy measurements

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Device Model : Park System XE Series

It is a scanning force microscope used for the optical, electrical or magnetic examination of the material surface, which provides high resolution, three-dimensional viewing of the surface with the help of a pointed tip of anatomic size.

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It performs applications such as evaluating the interactions between a particular atom and its neighboring atoms and examining the changes in physical properties resulting from atomic arrangements made by atom manipulation.

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100 micrometre length samples can be scanned.

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Maximum Z scan range can be extended from 12 μm to 25 μm .

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Device Model : ZEISS GeminiSEM 500

GeminiSEM 500, which includes InLens, EsB, EBSD, EDX, SE, BSD and STEM Detectors, has the feature of providing much better resolution even at low voltages.

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It is an electron microscope that provides high resolution images at low voltages with field emission gun. EDS and EBSD apparatus provides high precision surface characterization. It is based on the principle of taking images by converting electron-sample interactions into signals through detectors as a result of sending electron beam to the sample surface.

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Determination of crystallographic structure, phase analysis.

2

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EDS Detector

Kalitatif ve kantitatif elementel analiz

Qualitative and quantitative elemental analysis.

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Characterization of polymers, thin film and biological samples.

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Device Model : HITACHI SU1510

The Scanning Electron Microscope (SEM) providing micro or nano size observations of the surface of the specimen is imaging equipment that focuses accelerating electrons in a vacuum, and scans a sample surface with an electron beam, detects signals produced by the specimen at that time, and transforms the produced signals on the monitor screen.

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Electrons emitted from an electron source (Tungsten filament) focus on the sample surface through a condenser and objective lenses. Primary electrons interact with atoms on the surface of the material and produced secondary electrons that provide an examination of the microstructure morphology, and backscattered electrons that give compositional contrast, and characteristic X-ray signals that offer elemental analysis and then these signals are detected by detectors.

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Powder morphology, fracture surfaces and microstructure of metals, ceramics and composites, organic materials, porous structures, polymers

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Composite materials, powder metallurgy, surface characterization

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Surface coated materials, thin films, metal, ceramic, polymer, composite materials

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Device Model : Renishaw Invia Reflex

One-Point, depth, different resolution mapping and 3-dimensional (volume scanning) analyzes can be performed. Due to lasers with two different wavelengths, i.e. 532 nm and 785 nm, it is possible to analyze different types of materials effectively.

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Raman spectroscopy is based on measuring the inelastic scattering of light (Raman scattering). In this method, the differences (i.e. Raman shift) between the wavelength of the light interacting with the molecule and the wavelength of the light scattered are measured. If infrared (IR)-inactive molecules exhibit polarization change, then they are called as Raman-active. Therefore, infrared spectroscopy and Raman spectroscopy are complementary.

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Structure analysis in the development process of the battery, chip, etc

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Analysis of natural and synthetic materials

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Surface-Enhanced Raman Spectroscopy (SERS) in the development of various sensors and biosensors

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Device Model : Profilm 3D

It is used to analyze sample surface profiles to measure surface topography and film thickness for light reflective materials.

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The optical profilometer is a device that creates a three-dimensional surface profile for measuring the surface roughness of solid materials. It can also analyze the surface roughness of small areas or profiles. It can also measure the thickness of each layer separately in a few layers of transparently coated thin film.

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It has 3D surface scanning, roughness, topography and step height measurement options.

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White light interferometer (WLI) is used to measure surface profiles and roughness up to 0.05 micrometer. Its minimum vertical feature size dead-lines to 0.001 micrometer with PSI.

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Device Model : Emcotest Durascan G5

The resistance of the material against plastic deformation is determined using hardness measurement. Vicker and Brinell are the types of hardness in material science. While the Brinell hardness test can be mostly applied to metals, the Vickers test can be used for all metals and ceramics and has one of the largest scales among hardness tests. The term hardness may also refer to resistance to bending, scratching, abrasion or cutting.

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The hardness value is calculated depending on the size and depth of the notch on the sample surface. Hardness measurement consists of measuring the resistance of the material against the dipping of a tapered or spherical standard tip into the material. When the selected tip is immersed in the material under a load, it leaves a mark on the material. The hardness of the material is inversely proportional to this trace.

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Calculation of hardness value at different loads for different metals and their alloys.

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Determination of the hardness value of the phases in macro and micro scales.

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Determination of hardness values of high toughness materials.

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Device Model : Emcotest N3

The resistance of the material against plastic deformation is determined using hardness measurement. Performing the Rockwell test is generally easier and more accurate than other hardness test methods. The Rockwell test method is used on all metals, except when the test metal structure or surface conditions will bring too many variations; the recesses will be too large for the application; or restricts the use of sample size or sample shape.

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The Rockwell test measuring the depth of penetration of an indenter under a large load compared to the penetration made by a preload. The Rockwell scale is a hardness scale based on indentation hardness of a material. The Rockwell test measuring the depth of penetration of an indenter under a large load compared to the penetration made by a preload. There are different scales, denoted by a single letter, that use different loads or indenters.

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Metals showing spring effect are tested for hardness.

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Hardness determination of annealed metals

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Hardness determination of hard metals and steel sheets

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Device Model : Olympus / GX41

The GX41 inverted metallographic microscope is ideal for quick and reliable specimen assessment and for determining if metallurgical properties meet manufacturing specifications. In this microscope, which specializes in metal samples, phase, pore and grain boundaries and distributions can be viewed under a maximum of 100x magnification.

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With the help of light dropped on the metal sample surface, it is provided to display and size the sample surface by passing through the 5, 10, 20, 50 and 100 magnification lenses and 10 magnification eyepieces. Since the response of the phases and grain boundaries to the etching reagent is different, topography is created in contrast that will allow the surface to be displayed by etching.

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By examining the metal surface, evaluations are made for the metal sample.

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Along with the grain size and distribution of the sample, it is evaluated by assigning a scale according to certain specifications.

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Device Model : Shimadzu Servopulser

Fatigue testing is a specialised form of mechanical testing that is performed by applying cyclic loading to a bar or structure. These tests are used either to generate fatigue life and crack growth data, identify critical locations or demonstrate the safety of a structure that may be susceptible to fatigue

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In many applications, materials are subjected to vibrating or oscillating forces. The behavior of materials under such load conditions differs from the behavior under a static load. Because the material is subjected to repeated load cycles (fatigue) in actual use, designers are faced with predicting fatigue life, which is defined as the total number of cycles to failure under specified loading conditions.

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Determination of axial fatigue types that can occur in materials.

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Analysis of tensile, compression and three-point bending tests that determine the static behavior of materials.

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Device Model : Agilent 7890B

GC is a instrument used in the analysis of vaporizable and high temperature resistant samples. It is widely used for the diagnosis, determination and structure analysis of substances separated in the gas chromatography.

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In gas chromatography, after the substances in the mixture are separated from each other, ionization is determined depending on the time they leave the column.

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Raw material and quality control analysis

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Qualitative and quantitative analysis of petroleum components

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Qualitative and quantitative analysis in the field of Food and Environment (fatty acids)

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Device Model : Agilent 5977A MSD

GC / MS is a superior feature instrument used in the analysis of evaporable and high temperature resistant samples.

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It includes the steps of ionizing the argon plasma, which is brought to 10,000 K temperature by electromagnetic induction, separating the ionized elements by mass spectrometry and measuring the element concentrations by an electron multiplexing detector.

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1

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Chemical and Pharmaceutical Industry

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Qualitative and quantitative analysis of drug components.

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Qualitative and quantitative analysis of petroleum and oil components.

3

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Agriculture and Environment Applications

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Qualitative and quantitative analysis of pesticide etc. compounds

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Device Model : Agilent 7010 GC/MS Triple Quadrupole.

Agilent 7010B GC / MS Triple Quad is the most powerful trace detector that enables target analysis in complex matrices. MS/MS simplifies sample preparation and shortens the analysis process. In addition, it helps to provide high specificity and sensitivity for the analysis.

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It includes the steps of ionizing the argon plasma, which is brought to 10,000 K temperature by electromagnetic induction, separating the ionized elements by mass spectrometry, measuring the element concentrations by an electron multiplexing detector.

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Important advantages are qualitative and quantitative analysis.

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Qualitative and quantitative analysis of petroleum and oil components

3

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Agriculture and Environment Applications

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Device Model :

Agilent 1260 Infinity LC - Agilent 6460 Triple Quad

Both parent ions and product ions formed as a result of fragmentation can be traced by tandem mass spectrometry. Electrospray ionization (ESI) technique is used in the tandem mass spectrometer working with the Agilent 1260 Infinity LC system. Qualitative and quantitative analysis of the substances can be performed using LC-MS/MS even at very low concentrations.

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In chromatography, the components in the sample are separated from each other by moving at different speeds with the help of a mobile phase in a stationary phase. In mass spectrometry, these components are transformed into plus or minus charged ions, and the masses of these charged molecules or fragments of molecules are determined by the mass/charge ratio. While liquid chromatography separates the components in the sample, mass spectrometry allows the structure of each component to be determined.

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Qualitative and quantitative analysis of drugs, vitamins, etc.

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Qualitative and quantitative analysis of flavonoids, dyes, etc.

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Qualitative and quantitative analysis of amino acids, drug metabolites etc.

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Device Model : Thermoscientific Ultimate 3000

Compared to conventional HPLC, Ultimate 3000 UHPLC device has various advantages such as small size steel columns having smaller particle size, lower sample amount, higher separation power and shorter analysis time.

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It is a method based on the separation of the components dissolved in a liquid, leaving the column at different times by moving at different speeds by interacting with the stationary phase in the steel column under high pressure. It enables qualitative and quantitative analysis of the components.

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Qualitative and quantitative analysis of phenolic compounds, dyes, etc.

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Qualitative and quantitative analysis of proteins, enzymes, lipids, hormones etc.

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Device Model : Objet30 Prime

Small design and engineering teams can meet diverse material needs such as transparency, flexibility and biocompatibility in-house. The Objet30 Prime offers: 1. rigid materials in multiple opaque shades as well as clear, for beautiful detail visualization and prototypes that include see-through components; 2. Rubberlike materials for soft-touch features and flexible components; and 3. specialized materials such as High Temperature, Simulated Polypropylene and even Bio- compatible.

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Proven PolyJet 3D Printing is famous for smooth surfaces, fine precision and diverse material properties. It works a bit like inkjet document printing, but instead of jetting drops of ink onto paper, the print head jets microscopic layers of liquid photopolymer onto a build tray and instantly cures them with UV light. The fine layers build up to create a prototype or production part.

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28 microns (0.0011 in.) for Tango materials; 16 microns (0.0006 in.) for all other materials

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Device Model : Dimension 1200es Series

Dimension 1200es 3D Printers use ABSplusTM modeling material, a production-grade thermoplastic durable enough to perform virtually the same as production parts. Models printed with Dimension 3D Printers have customer-proven toughness — from commercial sprayers tested at pressures up to 60 psi, to final parts on M1 tanks normally machined in aircraft-grade aluminum.

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Dimension 1200es working with Fused Deposition Modeling (FDM) technology use the real ABSplus thermoplastic for production. Dimension 1200es 3D printer allows for good resolution or fast production options, with layer thicknesses of 0.254 mm (0.010 in.) or 0.33 mm (0.013 in.).

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Certificates

Dimension 1200es: CE/ETL
SCA-1200 destek temizleme
sistemi: CE

Dimension 1200es: CE/ETL
SCA-1200 support cleaning
system: CE

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Device Model : Artec Eva ve Artec Space Spider

The Artec Eva and Space Spider is the best choice for creating fast, textured and accurate 3D models of medium sized objects. It impresses with mobility, ergonomics and a slim housing. The Artec Eva and Space Spider delivers data at high speed, with precision and vibrant colors. The scanner eliminates the need for time-consuming calibration and the sticking of markers to the object at the beginning of each scan.

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A new and enhanced precision instrument for CAD users and engineers, Artec Space Spider is a high-resolution 3D scanner based on blue light technology. It is perfect for capturing small objects or intricate details of large industrial objects in high resolution, with steadfast accuracy and brilliant color.

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Up to 0.1 mm 3D resolution and up to 0.05 mm 3D accuracy, Wide range of applications in reverse-engineering, quality control and more

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Make fast, high-quality scans of any object and Hybrid geometry and color tracking and data capture

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Device Model : Beckman Coulter - Cytotflex B4R0V0

Flow cytometer is a cell analyzer that can measure from 4 fluorescent channels with 488 nm blue laser. The device can be upgraded to 4 laser systems. The flow cytometer system can analyze 6 different fluorescent parameters at once, including forward (FS) and side scatter (SS).

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Flow cytometry is a technique and apparatus based on a single pass along the flow channel and the basis for classification based on cell size and granularity in this order as a suspension of various cells. The surface and internal proteins of cells with this equipment, organelles and other components analysis and separation, size using laser and electronic technology are performed by the granularity and fluorescence emission as basis.

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The presence or absence of antigenic structures are determined by specific antibodies.

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It detects how much cells are present in which division phase of the cells.

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It identifies alive, dead and apoptotic cells.

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Cihaz Modeli : QuantStudioTM 3 Real-Time PCR System

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Device Model : QuantStudioTM 3 Real-Time PCR System

The Applied Biosystems® QuantStudio® 3 Real-Time PCR System is designed for users who need an affordable, easy-to-use real-time PCR system that doesn't compromise performance and quality. The simplified Design and Analysis software is ideal for all types of users.

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Real-time Polymerase Chain Reaction (PCR) is the ability to monitor the progress of the PCR as it occurs in real time. In real-time PCR (qPCR), reactions are characterized by the point in time during cycling when amplification of a target is first detected rather than the amount of target accumulated after a fixed number of cycles.

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PCR, Genotyping, Gene expression analysis, miRNA quantification, HRM, allele detection, copy number detection

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Pathogen detection and antibiotic resistance screening

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Transgene copy number detection, quantification of the residual DNA in final product therapeutics

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Device Model : BAB,BAB-TAM-F

It is a light microscope that allows fluorescent molecules (blue, red, green) to be displayed with inverted objectives with magnification of 4X, 10X, 40X, 40x (oil).

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The organelle, protein or molecule that is desired to be imaged in the translucent samples (cell, tissue, etc.) is treated with fluorescent dyes or antibodies, and after the fluorescence is excited in an inverted microscope, emission is displayed, the image is photographed and analyzed by camera.

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Investigation of fluorescent molecules in translucent samples, analysis of biological and chemical samples containing endogenous fluorescent molecules

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Cell viability analysis, visualization of protein expression and its analysis, determination of the movement of a molecule or drug in the cell by fluorescent staining.

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Determination of the effect of molecules on DNA, FISH studies, detection of mutations.

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Epoch™ 2 is a monochromator based micro volume spectrophotometer suitable for 6-384 well microplates, cuvettes and 2 µL measurements, and is a device that can work with as little as 0.5 µL. It is used in studies that require UV / Visible absorbance and photometric measurements, spectrum scanning, nucleic acid and protein concentration measurements, enzyme, apoptosis, GPCR and kinetic studies. Incubation, agitation and robot compatibility are standard features of the system.

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A microplate reader is a laboratory instrument that is used to measure chemical, biological or physical reactions, properties and analytes within the well of a microplate. A microplate consists of small wells in which separated reactions take place. These reactions convert the presence of an analyte or the progression of biochemical processes into optical signals. The microplate reader detects these signals and thus quantifies the parameter of interest.

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1

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*Clinical Medicine and Molecular
Biology*

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Used in ELISA, enzyme kinetics Cytotoxicity, Cell proliferation, total protein determination, nucleic acid purity analysis studies.

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Used in spectral scanning, food safety and quality, bacterial identification studies.

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Used in spectral scanning, drug production and drug discovery studies.

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Device Model : BioRad-ChemiDoc Imaging System

The ChemiDoc Imaging System provides fast, reliable, and sensitive imaging and documentation of electrophoretic gels and chemiluminescence western blots.

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This system is compatible with stain-free methodology, chemiluminescence detection, and a wide range of gel stains such as ethidium bromide, SYPRO Ruby, Coomassie, and silver stains.

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Identify and quantify proteins, nucleic acids, protein-ligand interactions, etc. with chemiluminescence, or colorimetric detection.

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proteins, including stain-free for total protein normalization, and nucleic acids (DNA/RNA) gel imaging and documentation.

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Qualitative and quantitative analysis of products in the biotechnological DNA, RNA and protein structure such as cloned gene products, drug delivery systems and therapeutic agents.

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Device Model : Thermo Scientific ASE 350

With the help of temperature and pressure, the analyte from the sample is extracted in a much shorter time and fully automatically using less solvent and sample compared to other extraction methods. In addition to providing high recovery, losses arising from individuals are also prevented.

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In the extraction method, the solvent and unwanted impurities are removed from the mixtures and the desired analyte molecule is separated from the solutions or solid mixtures. Extraction of organic compounds from various solid and semi-solid samples is provided automatically with temperature and pressure.

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Pesticides, polyaromatic hydrocarbons, polychlorinated biphenyls, dioxins and furans, total petroleum hydrocarbons.

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Pesticides, fats and lipids, acrylamide, antibiotic residues, antioxidants.

3

Kimya ve Farmasötik

Chemistry and Pharmaceuticals

Farmasötik ajanlar, flavonoidler, metabolitler.

Pharmaceutical agents, flavonoids, metabolites.

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Device Model : CEM Mars 6

Microwave systems are generally used to bring solid samples or inhomogeneous mixtures into solution before ICP-MS and ICP-OES measurements. The sample that is brought into solution with this system can be analyzed in ICP-MS or ICP-OES after this process.

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The transition of the sample to the solution phase in the microwave system is accomplished by the addition of strong acid or acid mixture, adjusting the pressure and temperature. Different sample preparation methods are applied according to the type of the sample and the elements to be analyzed.

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Fast and reliable measurements are provided by more than a hundred pre-programmed methods.

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Burning up to 40 samples with the same method at the same time.

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It is widely used in geology, environment, food, biology, chemistry etc.

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Device Model : OPT-PC-2C Plasma Cleaner

Plasma cleaning removes contamination on surfaces using energetic plasma or Dielectric Barrier Discharge (DBD) created from gas types.

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Argon and oxygen or mixtures such as air and hydrogen / nitrogen are used. Low pressure (around 1/1000 atmospheric pressure) plasma is generated using high frequency (> MHz to kHz) voltage to ionize gas.

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Instrument Properties and Application Area

1

Esneklik ve Hassas Kontrol

Flexibility and Sensitive Control

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Plasma cleaning systems are designed to remove nano-scale organic contamination on the wafer or substrate.

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Mixtures such as argon and oxygen or air and hydrogen / nitrogen can be used.

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For metallographic analysis, sample preparation steps must be etched with a suitable solution and after polishing of the cut at an appropriate speed in accordance with the cooling fluid such that the sample size is required for each sample. Molding bakelite process is more suitable for small samples polishing operations .

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In the sample preparation processes, firstly, large samples are cut in a closed cabin with the coolant in such a way that they do not change the phase. If the sample sizes are not suitable for obtaining a smooth surface or holding, they are molded by bakelite molding method. The sample surfaces are polished according to the roughness with suitable sanders and felts by rotating these sanders in the sample polishing device.

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For metallographic analysis, cutting the samples with suitable coolant at appropriate speeds

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Molding of the sample embedded in bakelite granules with temperature

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Grinding the sample surface with appropriate sands and felts according to the roughness by rotating these sands in the polishing device

A photograph of a modern building with a glass facade. The word "BITAM" is displayed in large, blue, illuminated letters on the upper part of the building. A semi-transparent blue rectangular overlay is positioned on the right side of the image, containing the word "Kurumsal" in white text. The building's glass reflects the sky and surrounding environment.

BITAM

Kurumsal

A photograph of a modern, multi-story building with a facade of large glass windows. The building is set against a sky with soft, colorful clouds, suggesting a sunset or sunrise. A semi-transparent blue rectangular overlay is positioned on the left side of the image, partially covering the building. The word "Institutional" is written in a white, italicized serif font across the middle of this blue overlay. In the foreground, several vehicles are parked along the street in front of the building, including a white sedan, a white hatchback, and two silver SUVs. A green tractor is visible on the right side of the building's ground floor.

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Important points regarding the samples sent to Necmettin Erbakan University Science and Technology Research and Application Center (BİTAM) for analysis are available at <https://www.erbakan.edu.tr/bitam/sayfa/7643>. Applications of samples that do not meet the sample acceptance criteria are not accepted by BİTAM. We kindly ask you to consider other important issues on the page.

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